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Sterile diaphragm seat valve

加工工艺 The processing technology

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The valve body adopts CNC precision machining, which ensures the anastomosis of the sealing surface of the valve cavity and the curve of the diaphragm, reduces the friction of the diaphragm and extends the use of the diaphragm.

The surface polishing of the valve cavity can be machined or electropolished according to the customer's requirements. The polishing degree can be reached $0.25\ \mu\text{m}$.

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Soft elastic material of diaphragm for those fibers, solid particles, the working medium of pollution such as catalytic material will not produce the sensitive reaction, in general also won't affect the valve working and sealing. Depending on the temperature of work or disinfection, and the chemical properties of the working medium, different materials can be selected.

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Diaphragm valve diaphragm is produced according to FDA standards.

We can also provide other standard certifications according to customers' requirements.

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Due to the different working conditions to choose different types of valves and material, before choosing body and diaphragm, must analyze the application of each product, especially for the use of chemical medicine, and high temperature caused by the chemical reaction. Through effective chemical data or expert certification, it can be tested to ensure the safety and long-term effectiveness of the product.

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
Typical diaphragm method for screw fixation, and opening a fixed form contrast, this way of fixed spread stressed area to bolt the whole on the surface, this is to prevent damage to the diaphragm vacuum conditions mechanical link.

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The aseptic diaphragm valve is recommended for $15\sim 30$ Angle (according to different specifications) according to its structural characteristics. It is helpful for the valve to be completely excreted after cleaning, not to cause the liquid to stay inside the valve.

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


The diaphragm is the most important part of the diaphragm valve

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The diaphragm is the only part of the body that is directly exposed to the working medium except for the valve body.

Diaphragm valves are separated from the working medium by the diaphragm by the actuator and the external air. In addition, diaphragm is also a moving part, and diaphragm valve is controlled by the upper and lower movement of the diaphragm to achieve the control and cutting of the working medium flow.

The KST diaphragm has been developed and tested for many years according to the test standard and the strict test specifications of the third party. These tests are as close as possible to the actual conditions of the site, and to conduct continuous tests under different specifications. For example, testing the diaphragm plate with an automated saturated vapor closed-loop system is one of them. The results of the test have a great influence on the shape design, the composition of the material, the body design, the actuator design, and the assembly of the valve.

In addition to the size of MA10 diaphragm is to rely on elastic rubber button directly connected with the valve actuator, the rest is all types of diaphragm by embedding the inside of a stainless steel studs connected to valve executive structure.

All sizes are of the same type of material diaphragm having the same actuator, and can be replaced with each other without changing the pressure block and stem of the valve.

Different sizes of valves have the same size of diaphragm, which can reduce the stock of spare parts.

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EPDM is a kind of diaphragm material specially developed. In the middle, the layer of vulcanized fiber enhances its mechanical strength and is always made into a membrane system. The structure of the diaphragm has high stability when temperature rise and pressure increase. In addition, the reinforced fiber layer is hardened on the embedded stainless steel bolt to ensure that elastic rubber and metal can be linked better. Therefore, EPDM diaphragm is more suitable for use in a vacuum.

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PTFE (PTFE) diaphragm has high anti-chemical corrosion, high stability, good elasticity, long life, high fatality, no cold deformation; And excellent performance not affected by temperature fluctuation, whether it is heat cycle, cold cycle, or steam sterilization system.

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The diaphragm of the MA10 specification is usually designed to be monolithic, and PTFE is glued to the back of EPDM.

The diaphragm is always made into a membrane system. This single-chip diaphragm has a small surface area, which is suitable for the short stroke valve and has excellent operation Canon, which is suitable for the high frequency of opening and closing frequency, such as filling.

MA10 specifications of the diaphragm by elastic rubber buttons and valve actuator connection, MA20 specifications of the diaphragm is through the embedded stainless steel studs is connected to the valve actuator, to eliminate the diaphragm single point in the center of the strong force, so as to prolong the service life of diaphragm piece.

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The diaphragm of MA25 to MA80 specification is two separate types, consisting of a single PTFE and a piece of EPDM. EPDM is supported on the back and the diaphragm is always made into a film closure. The diaphragm design is a bit is that when the valve is closed, the diaphragm is in the normal shape and body sealing, thereby significantly reducing the stress of the diaphragm, prolong the service life of the diaphragm.

For the two separate PTFE diaphragm, its stainless steel cylinder head bolts are embedded in the PTFE diaphragm. In order to eliminate the single point force of the diaphragm center, the flexible diaphragm suspension structure is used to connect the valve actuator.

MA	10-80	10-20	25-80
Material	EPDM	PTFE/EPDM	PTFE/EPDM
Design form	Single-chip membrane system open	Single-chip membrane system open	Two type membrane system close
Temperature	-40 ~ 150*	-20 ~ 150	-20 ~ 160

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The standard of valve body is 1.4435/316l stainless steel, which conforms to ASME BPE standard and EN10204 test certificate 3.1 requirement/material inspection report. All body parts are printed with a steel plate, which can be used to trace the characteristics of the body material and its physical composition.

The interior profile of the valve body is specially designed to meet the requirements of GMP specifications. The best self-cleaning capacity has no dead zone design. The remaining space of residual fluid was eliminated and the service life of the diaphragm was improved.

阀体分类 Classification of the body	接口尺寸 Interface size	制造工艺 Manufacturing process
两通阀体 Two-way valve body	10 ~ 80 mm / 3/8 ~ 3" 10 ~ 80 mm / 3/8 ~ 3"	锻造 Forging 铸造 Forging
多通道阀体 Multi-channel valve body	20 ~ 65 mm / 3/4 ~ 2"	整块钢机加工 The whole steel machining
罐底阀 Tank bottom valve	20 ~ 50 mm / 3/4 ~ 2"	铸造 Forging

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Casting stainless steel body is made of die casting, and with the shape of the final valve body, the wax mould can produce the same products.

By immersing the body wax film in the ceramic, the entire body wax film surface is covered with a solid ceramic.

Dissolve the body wax in the ceramic cavity, then inject the melted stainless steel solution. After cooling, knock off the coated ceramic, and the blank body will be finished.

The casting body size is very accurate and the outer surface is quite smooth and tidy.

In order to obtain high quality casting, the mold of the product is specially designed and handled optimally.

Considering the material structure and density problem, the valve body of KST is required by the detailed test data, which is strictly inspected to ensure the high reliability of the product quality.

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
The following table is: the diameter difference of the stainless steel pipe of different international standard sanitary grade with the example of DN25

ISO 1127	DIN 11850	DIN	ASTM269	BS O.D.	SMS	JIS G
(DIN 11866 系列 B) (DIN 11866 series B)	系列 1 Series 1	选择系列 Choose the series	ASTMBPE (DIN11866 系列 C)	4825	3008	3447
	系列 2 Series 2		ASTM269			
	系列 3 Series 3		ASTMBPE (DIN11866 系列 C)			
	(DIN11866 系列 A) (DIN11866 Series A)		ASTM269			

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
According to the international standards of several commonly used, the size of the valve diameter and wall thickness. These standards and sizes shown in the following table.

Install a plumbing system, in order to correctly choose consistent international pipeline interface standard is very important, so as to ensure the process piping system conform to the requirements of the health level.

If the valve interface standard and system pipeline and system interface standards are inconsistent, is easy to cause deformation within the interface or the steps, there is no guarantee that the emptying function of the valve.

Butt welding is the most common way of connection interface, without any other materials. Usually includes automatic welding and argon arc welding way, besides this way, we can also provide other way of connection.

h1= 铸造阀体 Forging body
h2= 锻造阀体 Forging body

尺寸 size : mm

对焊接口标准 Buttt welding interface standards					ISO 1127		DIN 11850		DIN 系列选择 DIN Series of choices	ASTM 269 ASME BPE	SMS 3008	JIS G 3447		
DN	NPS	MA	L(min)	L ₁	h ₁	h ₂	Φd×s	系列 1 Series 1	系列 2 Series 2	Φd×s	Φd	s	Φd×s	Φd×s
手动隔膜阀 气动隔膜阀 Manual diaphragm valve Pneumatic diaphragm valve														
10	3/8	10	25	72	9	9	17.2x1.6	12x1.0	13x1.5	-	-	-	-	-
15	1/2	20	25	120	13	16	21.3x1.6	18x1.0	19x1.5	-	12.7	1.65	-	-
20	3/4	20	25	120	16	16	26.9x1.6	22x1.0	23x1.5	-	19.05	1.65	-	-
25	1	25	25	120	19	19	33.7x2.0	28x1.0	29x1.5	-	25.4	1.65	25.0x1.2	25.4x1.2
32	1 1/4	40	25	153	24	26	42.4x2.0	34x1.0	35x1.5	-	31.75	1.65	33.7x1.2	31.8x1.2
40	1 1/2	40	25	153	24	26	48.3x2.0	40x1.0	41x1.5	-	38.1	1.65	38.0x1.2	38.1x1.2
50	2	50	30	173	32	32	60.3x2.0	52x1.0	53x1.5	-	50.8	1.65	51.0x1.2	50.8x1.5
65	2 1/2	50	30	173	32	32	-	-	-	-	63.5	1.65	63.5x1.6	63.5x2.0
65	2 1/2	80	25	216	47	47	76.1x2.0	-	70x2.0	-	63.5	1.65	63.5x1.6	63.5x2.0
80	3	80	30	254	47	47	88.9x2.3	-	85x2.0	-	76.2	1.65	76.1x1.6	76.3x2.0

MA 表示膜片大小 MA said the size of the diaphragm

根据要求，也可以提供其他的对接标准 According to the requirements, can also provide other standard docking

上表中黑体所示的尺寸为首选规格 The above table shown in bold in the size of the preferred specification

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
Clamp interface is one of the most common way of connection, the valve in the pipeline easily installed and remove. Clamp interface is generally designed to connect face to face, to prevent leakage, no cracks.

Clamp face has a machined into the annular groove, when used with a special molding of EPDM or PTFE gasket to seal material.

Washer is placed in the middle of two clamp, use a chuck clamping band, make whole interface is completely sealed.

Generally speaking, the clamp on the butt welded to the valve interface, then according to the requirements of the valve body surface polishing and polishing.

After welding, 100% of the hoop end passes through the visual inspection and pressure test, for all the international standard pipe, can be welding corresponding card hoop.

If the connection card hoop don't match each other, may lead to deformation of interface or the steps, there is no guarantee that the emptying function of the valve.

If it can be properly installed with a clamp valve, after installation of piping systems will be smooth, no cracks, and self aligned connection, thus reducing the pipeline was the risk of contamination, maximum limit reduces the medium disorder and system pressure drop.



卡箍接口标准 管径接口标准 The clamp mouth interface standards Pipe diameter interface standards			类似于 ISO 2852 ISO 1127			DIN 32676 DIN 11850			ASME BPE ASME BPE			ASME BPE ADME BPE			SMS 3017 SMS 3008		
标准 FIF			DIN EN 558-1			DIN EN 558-1			DIN EN 558-1			ASME BPE DT-V-1			DIN EN 558-1		
DN	NPS	MA	L3	d2	d1	L3	d2	d1	L3	d2	d1	L3	d2	d1	L3	d2	d1
8	1/4	10	*63.5	10.3	25.2	-	-	-	-	-	-	63.5	4.57	25.2	-	-	-
10	3/8	10	-	-	-	*63.5	10.0	34.0	-	-	-	63.5	7.75	25.2	-	-	-
15	1/2	20	108	18.1	50.5	108.0	16.0	34.0	108.0	9.40	25.2	101.6	9.40	25.2	-	-	-
20	3/4	20	117	23.7	50.5	117.0	20.0	34.0	117.0	15.75	25.2	101.6	15.75	25.2	-	-	-
25	1	25	127	29.7	50.5	127.0	26.0	50.5	127.0	22.10	50.5	114.3	22.10	50.5	127.0	22.6	50.5
32	1 1/4	40	146	38.4	50.5	146.0	32.0	50.5	-	-	-	-	-	-	146.0	31.3	50.5
40	1 1/2	40	159	44.3	64.0	159.0	38.0	50.5	159.0	34.80	50.5	139.7	34.80	50.5	159.0	35.6	50.5
50	2	50	190	56.3	77.5	190.0	50.0	64.0	190.0	47.50	64.0	158.8	47.50	64.0	190.0	48.6	64.0
65	2 1/2	80	216	72.1	91.0	216.0	66.0	91.0	216.0	60.20	77.5	*222.3	60.20	77.5	216.0	60.3	77.5
80	3	80	254	84.3	106.0	254.0	81.0	106.0	254.0	72.90	91.0	222.3	72.90	91.0	254.0	72.9	91.0

手动隔膜阀 DN10-20mm(3/8"-3/4")
Manual diaphragm valve DN10-20mm(3/8"-3/4")

产品特点 Product features

- 升降式手轮
- 带就地阀位指示的密封阀盖
- 可调节的闭合限位
- 阀体隔膜和阀体间具有固定的密封圆环
- 柔性隔膜悬浮结构
- Lift handwheel
- Seal valve cover with local valve position indication
- Adjustable closing limit
- The body diaphragm and valve body have a fixed seal ring
- Flexible diaphragm suspension structure

控制方式: Control mode	手动 Manual
最大工作压力: Maximum operating pressure	10 bar
最高工作温度: Maximum operating temperature	160°C
隔膜材质: Diaphragm material	EPDM、PTFE
阀体材质: Valve body material	锻造不锈钢 Wrought stainless steel 1.4435 / 316L、ASME / BPE 铸造不锈钢 Wrought stainless steel 1.4435 / 316L 其它合金 Other alloy
连接方式: Connection method	对焊接口 Butt welding interface 卡箍和法兰接口 Clamp and flange interface 特殊接口, 请咨询 Special interface please consult
阀盖适用于: Cover is suitable for	两通阀体 Two-way valve body 焊接组合式阀体 Welded sectional body T-型阀体 T type valve body 多通道阀体 Multi-channel valve body 罐底阀体 Tank bottom valve body

手动隔膜阀 DN25-80mm(3/8")
Manual diaphragm valve DN25-80mm(3/8")


产品特点 Product features

- 不锈钢阀盖和耐高温工程塑料手轮
- 带就地阀位指示的止升式手轮
- DN50 以下的阀体隔膜和阀体间具有固定的密封圆环
- 柔性隔膜悬浮结构
- 包覆式安装的阀体隔膜
- Stainless steel valve cover and high temperature resistant engineering plastic hand wheel
- With valve position indicator on the spot check up the handwheel
- (the following valve body between the diaphragm and the valve body has a fixed seal ring
- The flexible diaphragm suspended structure
- The diaphragm body of cladding type installation

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
可选附件 Optional accessories

- 可选装多种控制附件，形式多样
- 进气口方向可与工作介质流向成 90°
- 高压消毒
- Optional multiple control accessories, various forms
- Air inlet direction to the flow direction of medium can work with 90°
- High pressure disinfection

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Pneumatic diaphragm valve DN25-80mm(1"-3")

产品特点 Product features

- 高循环的活塞式不锈钢气动执行机构
- 设计紧凑，气动头的外径尺寸与连接隔膜的阀盖法兰尺寸一样
- 非常适合于多通道阀体和多阀门的组合安装
- 进气口方向与工作介质流向一致
- 阀体隔膜和阀体之间具有固定的密封圆环
- 柔性隔膜悬浮结构
- 包覆式安装的阀体隔膜
- 清洁、抛光的执行机构外表面非常适合消毒后冲洗
- High cycle piston stainless steel pneumatic actuator
- The design is compact, and the outer diameter of the pneumatic head is the same as the valve cover flange of the diaphragm
- Suitable for multi-channel body and multi-valve assembly installation
- The inlet direction is in line with the working medium
- A fixed seal ring between the valve body diaphragm and valve body
- Flexible diaphragm suspension structure
- Covered body diaphragm
- The outer surface of the clean and polished actuator is well suited for post-sterilization

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
Different actuators can achieve higher working pressure

规格 (DN)	外形尺寸 (mm)							
	MA	L	L ₁	AxB	H ₁	H ₂	H ₃	D
25	25	25	120	73x79	146	66	133	75
40	40	25	153	96x105	180	75	160	105
50	50	30	173	111x130	216	77	180	105
65	80	30	216	190x170	309	135	285	175
80	80	30	254	190x170	309	135	285	175

气动隔膜阀 DN25-50mm(1"-2") Pneumatic diaphragm valve DN25-50mm(1"-2")

产品特点 Product features


- 耐高温工程塑料的活塞式气动执行机构
- 设计紧凑
- 工程塑料执行机构具有高抗热传递
- 进气口方向与工作介质流向一致
- 阀体隔膜和阀体之间具有固定的密封圆环
- 柔性隔膜悬浮结构
- 包覆式安装的阀体隔膜
- 清洁光滑的执行机构外表面非常适合冲洗
- Piston pneumatic actuator for high temperature engineering plastics
- Compact design
- Engineering plastic actuator has high heat resistance transmission
- The inlet direction is in line with the working medium
- A fixed seal ring between the valve body diaphragm and valve body
- Flexible diaphragm suspension structure
- Covered body diaphragm
- The clean, smooth surface of the actuator is well suited for flushing

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
Different actuators can achieve higher working pressure

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Pneumatic diaphragm valve DN08-50mm

型号	d	D	D1	L	A	H	执行器
DN08	8	25.2	20.2	89	64	122	050
DN10	10	25.2	20.2	89	64	122	050
DN15	15	43.5	50.5	108	80	145	063
DN20	19	43.5	50.5	118	100	167	080
DN25	22.6	43.5	50.5	130	100	170	080
DN32	28	43.5	50.5	150	100	175	080
DN40	36	43.5	50.5	161	126	225	100
DN50	50	56.5	64	172	126	238	100

产品特点 Product features

- 高循环的活塞式不锈钢气动执行机构
- 设计紧凑，气动头的外径尺寸与连接隔膜的阀盖法兰尺寸一样
- 非常适合于多通道阀体和多阀门的组合安装
- 进气口方向与工作介质流向一致
- 阀体隔膜和阀体之间具有固定的密封圆环
- 柔性隔膜悬浮结构
- 包覆式安装的阀体隔膜
- 清洁、抛光的执行机构外表面非常适合消毒后冲洗
- High cycle piston stainless steel pneumatic actuator
- The design is compact, and the outer diameter of the pneumatic head is the same as the valve cover flange of the diaphragm
- Suitable for multi-channel body and multi-valve assembly installation
- The inlet direction is in line with the working medium
- A fixed seal ring between the valve body diaphragm and valve body
- Flexible diaphragm suspension structure
- Covered body diaphragm
- The outer surface of the clean and polished actuator is well suited for post-sterilization

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
T - type valve body is welded to a circular tube in a liquid or pure water system. This technique can be used for sampling or water use without a dead Angle. In addition, according to the requirements of the process, the mixing of two media can also be carried out. It is important that the seal of the t-type valve should be as close to the main pipe as possible.

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
**T 型阀或 U 型阀 (零肆角阀)
T type valve or u-type valve (zero Angle valve)**

1 个用水点阀座
1 water point seat

安装建议:
S3 口朝下
Installation Suggestions:
S3 mouth down

规格		外形尺寸 (mm)								
DN1	DN2	L	L1	L2	L3	H	H1	H2	B	B1
20	20	93	16.7	25	59	36	22.3	19.0	124	25
25	20	96	20.1	25	59	42	25.3	19.0	124	25
25	25	96	20.1	25	59	42	25.3	19.0	124	25
40	20	108	26.7	25	59	54	31.3	19.0	124	25
40	25	108	26.7	25	59	54	31.3	19.0	124	25
40	40	122	28.8	25	71	57	34.2	26.0	152	25
50	20	121	33.4	25	59	66	37.3	19.0	134	30
50	25	121	33.4	25	59	66	37.3	19.0	134	30
50	40	135	35.5	25	71	69	40.2	26.0	162	30
50	50	151	36.7	30	85	71	42.3	32.0	184	30

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
The tank bottom valve is mainly used in biological pharmacy, food and beverage and cosmetics, etc. The inner surface is smooth and free of dead zones, preventing the retention of the medium and avoiding the potential process pollution. The tank bottom valve has the same characteristics and functions as the standard diaphragm valve and USES the same valve components, such as actuator, diaphragm, etc. Tank bottom valves are usually soldered directly to the bottom of the tank and can be soldered to the tank wall as a tank wall valve or sample valve.

The typical application of tank bottom valve is to realize the discharge, emptying, sampling, cleaning, disinfection, washing and separation of the material in the tank.

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The standard design of the bottom valve is to carry only one seat exit. In addition, many functional forms of tank bottom valves are available for different applications, such as sampling, disinfection, and multiple export structures to standardize these products.

The preferred way to install the bottom valve is to weld it directly to the jar. The direct installation of the valve and tank greatly reduces the amount of the medium, which is the most important standard for using a tank bottom valve. However, if the tank bottom valve needs to be removed from the tank, a flange connection or clamp connection may be selected between the valve and the tank.

The interface of the valve outlet of the bottom valve is: welding connection, health card coupling.

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Intelligent valve locator is designed for one-piece pneumatic valve, especially for Angle seat valve and diaphragm valve products. Its products can be divided into process control valve and proportional control valve.

This product is easy to operate and rich in software functions. It can be operated easily through LCD panel and key. The product senses valve opening through position sensor and minimizes energy consumption.

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The design of the welded sectional diaphragm valve is based on GMP specification to reduce the dead Angle and improve the process flow of sanitary production equipment. Welding combined diaphragm valve can be made very simple single valve and pipe combination, can also make it very complicated, different size of the valve body welding together trip valve group of various valves. All forms of welding interfaces can be made into solder composite diaphragm valves.

Each welded composite diaphragm valve has strict quality control. Each weld on the valve is polished according to the inside surface of the valve to ensure the internal surface quality of the valve.

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
D - the rule is a rule that reflects the dead Angle based on the relationship between size B and D described by ASME BPE standard. The rule defines a very useful principle that describes the maximum dead Angle that is allowed between interconnected parts of the plumbing system. The dead corner stone is reflected by the absolute value of B or the ratio of B/D.

The ratio of B/D should be converted between 2:1 and 5:1 depending on the combination of the size of the body and the size of the valve body. If the value has been specified and the solder combined diaphragm valve cannot meet this requirement, the multi-channel valve can be used to solve this problem.

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1) GMP valve

GMP valve (abbreviation for Good Manufacturing Practice) is mainly used in high purity water or other distribution system with water point vertical pipe, which can effectively reduce the dead Angle of pipeline. There are two types of valves designed by GMP: one is the zero Angle valve with a 90-degree Angle bend, and the other is the zero Angle valve connected to the valve. In the valve phase connection, one of the valves is placed horizontally and maintains a certain self-discharge Angle, and the other does not have a bacterial breeding or process cross infection sampling port.

For this type of structure, whether the main valve or the L valve or the outlet, the maximum size of the unit is DN80 (3 ').

Here are some possible GMP connections:

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2) SAP valve

SAP valve (the abbreviation of Sterile Access Port) is mainly used in horizontal piping system, including the main valve is in horizontal, and maintain the emptying of Angle, via hole in the main valve of the lowest waterway discharge point. The SAP structure can be a combination of main and partial valves, as well as a combination of valve and vertical or horizontal valves.

For this kind of structure, whether the main valve or the channel valve or the branch, its maximum size can reach DN80 (3 ').

Here are some possible SAP connectivity structures:

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According to the requirement, the pneumatic head can install the mechanical limit and the signal feedback device simultaneously.